


OMNI-STAT®

CELOX™ Rapid


How to Use/Recommendations for Removal



1a Identify the point of bleeding and apply direct pressure.



1b Swab excess blood where practical.




2 Peel open pouch. Remove the Celox Rapid and take one end of Z-folded gauze.




3 Tightly pack the unfolding Celox Rapid directly to the source of the bleeding, feeding in the Celox Rapid until the wound is packed above skin level. Excess Celox Rapid can be cut or torn if necessary.



4 Apply FIRM pressure directly to the wound for 60 seconds. If bleeding persists, apply additional pressure until bleeding stops.



5 Secure with a dressing to maintain pressure on the wound.



6 Discard any remaining Celox Rapid.

REMOVAL INSTRUCTIONS:

At first dressing change after use, hemostatic gauze should be cleansed and removed from the wound using standard wound cleansing protocols. If required, should be soaked with saline prior to removal (physically) and then any residual irrigated away with water or saline.

DURATION OF USE:

The device and residuals should be removed from the wound within 72 hours from application.

For full details, warnings and removal instructions refer to manufacturer's instructions for use.

1. Kauvar, D.S., Lefering, R., and Wade, C.E. (2006) Impact of hemorrhage on trauma outcome: an overview of epidemiology, clinical presentations, and therapeutic considerations. *J Trauma* 60, S3-11.
Copyright © 2024 Medtrude Products Ltd. CELOX™ is a trademark of Medtrude Products Ltd. OMNI-STAT® is a registered trademark of Omni-stat Medical Inc. All rights reserved.

PRODUCT EVALUATION FORM

Evaluation forms can be emailed to info@omni-stat.com or given to your local sales representative

Facility/Hospital Name: _____ Date: _____

Provider Name: _____ Provider Signature: _____

Product Preparation Used

- OMNI-STAT Granules (3g) CELOX™ Z Fold Gauze (5ft) CELOX™ Granules (15g)
 OMNI-STAT Gauze (4in x 4in) CELOX™-A (6g) CELOX™ RAPID

Type of Wound

- Gunshot
 Laceration/Abrasion
 Puncture
 Blunt Trauma
 Skin Tear
 Diabetic Foot Ulcer
 Pressure Ulcer
 Other: _____

Hold Pressure Time

- 1 minute
 3 minutes
 5 minutes
 10 minutes
 Other: _____


Was the Patient Anticoagulated?

- YES,
If Yes, please list type
of anticoagulant: _____
 NO

Was hemostasis achieved?

- YES
 NO

Wound Location



Please circle the location of the wound on the diagram above

Observations/Results

